

## Analog Discovery 2 Complete Control Panel – User Guide

This application provides a **graphical interface** for controlling and monitoring the [Diligent Analog Discovery 2](#) device. It allows you to use the device's oscilloscope, function generator, power supplies, digital I/O, data logger, spectrum analyzer, protocol analyzer, and network analyzer features.

---

### Getting Started

#### Requirements

- **Hardware:** Diligent Analog Discovery 2
- **Software:**
  - Windows OS
  - [Diligent WaveForms SDK](#) (dwf.dll must be available)
  - Python 3.x
  - Required Python packages: [tkinter](#), [numpy](#), [matplotlib](#)

#### Launching the App

1. Connect your Analog Discovery 2 to your PC via USB.
2. Ensure [dwf.dll](#) is accessible (in the same directory or in your PATH).
3. Run the app:

```
python complete_ad2_gui.py
```

- - 
  - 
  -
4. The main window will appear.
- 

### Main Features & Tabs

#### 1. Oscilloscope & Function Gen

- **Connect/Disconnect:** Use the buttons at the top to connect or disconnect from the device.
- **Oscilloscope:**
  - Enable/disable channels, set voltage range, timebase, and trigger.

- Start/Stop/Single acquisition.
    - Save acquired data as CSV.
    - View real-time plots of both channels.
  - **Function Generator:**
    - Enable/disable each channel.
    - Select waveform (Sine, Square, Triangle, DC, Sawtooth, Noise).
    - Set frequency, amplitude, and offset.
- 

## **2. Power Supply**

- **Enable/disable positive and negative supplies.**
  - **Set voltage for each supply.**
  - **Monitor current draw in real time.**
  - **Enable/disable all supplies at once.**
- 

## **3. Data Logger**

- **Select channels to log.**
  - **Set logging interval and duration.**
  - **Choose output CSV file.**
  - **Start/Stop logging.**
  - **Monitor progress and status.**
- 

## **4. Spectrum Analyzer**

- **Set start/stop frequency and number of samples.**
  - **Start/Stop spectrum analysis.**
  - **View real-time spectrum plot (simulated data).**
- 

## **5. Digital I/O**

- **Set digital output states for all 16 DIO pins.**
- **Monitor digital input states in real time.**

---

## 6. Protocol Analyzer

- **Select protocol (SPI, I2C, UART, CAN).**
- **Set clock and data pins.**
- **Start/Stop protocol analysis (feature not implemented, placeholder only).**
- **View decoded data (not implemented).**

---

## 7. Network Analyzer

- **Set frequency sweep parameters (start/stop freq, points, amplitude).**
- **Start/Stop frequency sweep.**
- **View Bode plot (magnitude and phase) of S11 response.**
- **Export data to CSV or save plot as image.**

---

## 8. Settings

- **View device information.**
- **Calibrate Oscilloscope and Function Generator.**
- **Reset device to factory defaults.**
- **Set advanced parameters (buffer size, timeout).**

---

## General Usage Tips

- **Connect the device** before using any features. Most controls are disabled until a device is connected.
- **Calibration and reset** operations are available in the Settings tab.
- **Data export** is available for Data Logger and Network Analyzer.
- **Some advanced features (Protocol Analyzer, real hardware spectrum/network analysis)** are placeholders or simulated. You may need to implement or extend these for full hardware support.

---

## Troubleshooting

- **WaveForms library not loaded:** Ensure [dwf.dll](#) is present and accessible.

- **Device not found:** Check USB connection and drivers.
  - **Feature not implemented:** Some tabs (e.g., Protocol Analyzer) are placeholders and will show an info dialog.
- 

### Closing the App

- Simply close the window or use the standard OS controls.
- 

### Support

For hardware issues, refer to the [Diligent Analog Discovery 2 documentation](#).  
For software issues, check your Python environment and dependencies.

---

**Enjoy exploring your Analog Discovery 2!**